



Gold Hpdry

Rapid urease test for *Helicobacter pylori* detection

- ▶ Detection of *Helicobacter pylori* in gastric mucosal bio-samples.
- ▶ Easy, fast and reliable, it is performed in the endoscopy suite and allows rapid diagnosis of *Helicobacter pylori* infections.



Test in 4 steps:



Reference	Description	QTY
94000000	Rapid urease test for Helicobacter pylori detection	50

Characteristics

- ▶ **Rapid detection** of Helicobacter pylori infection:
 - First reading at **5 minutes**.
 - Final reading at **60 minutes**.
- ▶ **Variation in colour intensity** to assess the **degree of infection**.
- ▶ **98% sensitivity** and **97% specificity**.
- ▶ Can be stored for **2 years without any particular storage method**.
- ▶ **Economic**.

Description: Urease rapid test for the detection of Helicobacter pylori GOLD Hpdry. **Indications:** Detection of Helicobacter pylori in gastric mucosal biopsies, based on the detection of urease activity. **Class:** I. **Notifying body:** Self-certification. **Manufacturer:** GrupoBios S.A. **Proper use:** Before use, carefully read the instructions for use or the labelling of the device. **Version :** GHP_EXP_2 Rev. 02 Date 22-09-22.

Creo Medical authorised local distributor

Creo Medical S.L.

Pol. Ind. Cordovilla, Calle D. N°1
31191 Cordovilla, Navarra (España)
Tel.: +34 948 293 324
E-mail: info.es@creomedical.com

