


ScopeGuard


Biteblocks



Large central opening
Allows the passage of all types of endoscopes



Small size model
For patients with limited mandibular mobility and Infants




Internal dental stops
Avoiding the expulsion of the mouthguard by the patient

Lateral indentations
Facilitate the guiding of the endoscope during its introduction

The most complete range of Biteblocks suitable for all upper digestive endoscopies



Models with silicone dental protection
Protecting the patient's teeth



Oxygenating biteblock
For endoscopies where O₂ is required

Eco model
Taking into account budgetary constraints





Single Use - Individually Packed supplied with Tracability Labels

Standard Biteblocks - Non Sterile

	Reference	Description	Central opening (mm)	Qty
	AN10018000	ScopeGuard Universal Biteblock	27 x 22	100
	AN2360000	ScopeGuard Plus Eco Biteblock + Head strap	26 x 20	100
	AN10018002	ScopeGuard Plus Universal Biteblock + Head strap	27 x 22	100
	AN10018003	ScopeGuard Plus Mini Small size Biteblock + Head strap	20 x 16	100
	ENDO0209	ScopeGuard Plus Protect Biteblock with dental protection + Head strap	26 x 19	100
	ENDO0209A	ScopeGuard Plus Biteblock + Head strap	26 x 19	100
	AN10018011	ScopeFix Head strap	-	100

BiteBlocks oxygenators - Non Steriles

	Reference	Description	Central opening (mm)	Qty
	ENDO044	ScopeGuard O ₂ Biteblock for oxygen insufflation + Head strap	23 x 22	100
	ENDO044F	ScopeGuard O ₂ Protect Biteblock for oxygen insufflation with dental protection + Head strap	23 x 22	100

SG_E_4, Rev. 04 Date 16-09-22

Creo Medical authorised local distributor

Creo Medical UK Ltd

Chestnut House, N.Minster Business Pk.

YO26 6QR York, United Kingdom

Email: info.uk@creomedical.com

Tel: +44 (0) 1904 786888

Technical changes and mistakes reserved. Images may show non standard options.

