

UroAccess

Ureteral Access Sheath

4 configurations where the sheath is declined:

- ▶ In **2 diameters**: 10 - 12Fr and 12 - 14Fr.
- ▶ In **2 lengths**: 35 cm and 45 cm.



**Designed to create and maintain
secure access to the kidneys.**



UroAccess

Ureteral Access Sheath

Reference	Ø Sheath exterior (Fr)	Ø Sheath interior (Fr)	Length (cm)	QTY
UAS-101235-041	10	12	35	4
UAS-101245-042	10	12	45	4
UAS-121435-038	12	14	35	4
UAS-121445-026	12	14	45	4

Characteristics

- ▶ **Optimised shape of the dilator tip** for easy insertion.
- ▶ **Bevelled access sheath end** for atraumatic transition between dilator and sheath



- ▶ Hydrophilic coating provides **excellent glide** during insertion.
- ▶ **Reinforced** and **radiopaque** sheath.
- ▶ Sheath is **graduated** along its length.

- ▶ **Luer-Lock connection** for contrast medium injection.
- ▶ **Locking system** of the expander in the sheath for a reliable coupling.



Description: UroAccess urethral access sheath. **Indications:** Provides a channel for an interventional device to reach the area where the operation takes place. **Class:** IIa. **Notified body:** CE1984. **Manufacturer:** Plasti-Med Plastik Medikal Ürünler San.ve Tic.Ltd.Şti. **Proper use:** Before use, carefully read the instructions in the package insert or device labelling. **Version:** UA_E_1.Rev.01 Date 02-12-2021

Distributed in UK by:

Creo Medical UK Ltd.

Chestnut House, Northminster Business Park

YO26 6QR York, United Kingdom

Email: info.uk@creomedical.com

Tél: +44 (0) 1904 786888

