


WiScope[®]

Single-Use Digital Flexible Cystoscope

A high-performance device, always operational and available, that can be used in the operating suite, in consultation or at the patient's bedside



WiScope®


Single-Use Digital Flexible Cystoscope

Reference	Deflection	Insertion sheath Ø (Fr / mm)	Operating channel Ø (Fr / mm)	Sheath length (mm)	Compatibility		PKG	MOD
					Tablet	Processor		
00060011	Reverse	15,3 / 5,1	6,6 / 2,2	380	✓	✗	2	6
00060008	WiScope® Medical Tablet						1	1

Cystoscope with standard deflection also available, please contact us.

Characteristics

- ▶ Get **high-performance equipment** that is **always available and 100% functional** for your cystoscopic examinations, without delaying patient care. In addition to its immediate availability, the **WiScope®** single-use cystoscope:
 - Requires **no investment, no repairs, no maintenance and no microbiological control**,
 - Eliminates **decontamination and sterilisation** procedures,
 - Eliminates the **risk of cross-contamination**,
 - Delivers the expected performance in terms of **handling and vision quality**,
 - **Enables rapid intervention and emergency management.**
- ▶ The **WiScope®** medical tablet is **lightweight, portable and compact**, so it can be **easily transported from one department to another**, enabling examinations to be carried out at the patient's bedside.



Description: WiScope® Single-Use Digital Flexible Cystoscope. **Indications:** access, visualize, and perform procedures in the urinary tract and kidney. **Class:** IIa. **Notified body:** CE0123. **Product use:** Before use, read attentively the instructions on the Instructions for Use or on the product label. **Version :** WSC_E_3. Rev.03 Date 12-12-23

Creo Medical authorised local distributor:

Creo Medical UK Ltd

Creo House, Unit 2 Beaufort Park
Beaufort Park Way, Chepstow, NP16 5UH
Email : info.uk@creomedical.com
Tél : +44 (0) 1904 786888

