

ClearCap Distal Attachment



Reference	Applicable endoscope	O.D	I.D	Drainage holes	Number of Use	Validity of Sterilization
FM-EC0001	Lower GI	Ø 15mm	Ø 12.6mm	2	Single Use	3 years
FM-EC0002	Upper GI	Ø 11.8mm	Ø 9.4mm	2	Single Use	3 years
FM-EC0003	Upper GI	Ø 11.35mm	Ø 8.9mm	2	Single Use	3 years
FM-EC0004	Upper GI	Ø 12.4mm	Ø 10mm	2	Single Use	3 years
FM-EC0005	Lower GI	Ø 13.4mm	Ø 11mm	2	Single Use	3 years
FM-EC0006	Lower GI	Ø 14mm	Ø 11.6mm	2	Single Use	3 years
FM-EC0007	Lower GI	Ø 15.7mm	Ø 13.3mm	2	Single Use	3 years

Characteristics

- ▶ **Maintaining a Clear View:**
 - The ClearCap Distal Attachment helps maintain a clear view with colourless transparent material.
- ▶ **Protecting the tissues:**
 - Soft material facilitates endoscope insertability and reduces tissue trauma.
- ▶ **Two draining holes:**
 - Rapid drainage of fluid to improve procedural visibility.

Creo Medical authorised local distributor

Creo Medical UK Ltd

Chestnut House, N.Minster Business Pk.

YO26 6QR York, United Kingdom

Email: info.uk@creomedical.com

Tel: +44 (0) 1904 786888

