

ScopeInject

Injection Needles

Key Characteristics

Intended for endoscopic injection into GI mucosa

Scope Inject:

Ergonomic handle: designed to be operated with one hand

Designed for safety: keeps the needle in "extended" and "retracted" position

Transparent sheath to visualise any movement of the needle

Additional features of ScopeInject Plus:

Coated: intended to cause less friction and to prevent kinking

Metal end: designed for easy deployment through the channel

Designed for usability with configurations best suited for the clinical need

Sharp needle, strong puncture force

Variable gauge and needle lengths

Sterile and single use



ScopeInject					
Reference	Catheter		Needle		Box QTY
	Ø (mm)	Length (mm)	Ø Gauge	Length (mm)	
IN-241621045081	2.4	1600	21	4	10
IN-241623045082	2.4	1600	23	4	10
IN-241625045083	2.4	1600	25	4	10
IN-242321045081	2.4	2300	21	4	10
IN-242323045082	2.4	2300	23	4	10
IN-242325045083	2.4	2300	25	4	10

ScopeInject Plus					
Reference	Catheter		Needle		Box QTY
	Ø (mm)	Length (mm)	Ø Gauge	Length (mm)	
IN-2416-2304M-5313	2.4	1600	23	4	10
IN-2416-2504M-5313	2.4	1600	25	4	10
IN-2423-2304M-5313	2.4	2300	23	4	10
IN-2423-2504M-5313	2.4	2300	25	4	10

Discover more:



Creo Medical Ltd.

Unit 2, Creo House
Beaufort Park Way
Chepstow NP16 5UH, UK
+44 (0) 1291 637 300
customerservice@creomedical.com

Creo Medical Inc.

100 Reserve Road
Suite B400
Danbury, CT 06810, USA
+1 866-226-1170
UScustomerservice@creomedical.com

Creo Medical Pte Ltd.

60 Paya Lebar Road, #09
01/02/03/04 Paya Lebar Square
Singapore
409051

in [showcase/creo-endoscopy](#)

in [company/creo-medical](#)

X [@CreoMedical](#)

X [@CreoEndoscopy](#)

Technical changes and mistakes reserved. Images may show non standard options.

Creo Medical, the CM logo, Kamaptive and Speedboat are trademarks of Creo Medical Limited. Certain names and marks are registered in the UK and other territories. All cited trademarks are the property of their respective owners. © 2024 Creo Medical. All rights reserved.

MKT-GEN-003-004 | 2024-April

